

Certificate in relation to air transportation by passengers with walking disabilities

We hereby certify that the our product

smoov (one) O10

including the lithium-ion batteries meets the requirements of

UN Manual of Tests and Criteria" Part III subsection 38.3

and therefore are suitable for transportation by passengers according

IATA DGR Chapter 2.3.2.4

The rated energy of the drive unit battery **216 Wh** is below the IATA limit of 300 Wh. Also the battery of the control unit with 9,36 Wh is below the IATA threshold.

Since drive unit and control unit are portable and contain batteries, we recommend that you take them with you to the passenger cabin. Both units have to be switched off before boarding and protected from damage and unintentional activation.

We strongly recommend to contact the airline sufficiently in advance to request approval for transport.

This certificate does not apply to damaged or defective devices or batteries.

Due to the frequent amendments of the respective regulations, which partly occur on an annual basis this certificate is valid until **31st December 2020**.

Albstadt, 22nd of November 2019

Alber GmbH



Hansjörg Reiner
Managing Director