

Data sheet of the via GO V24 lithium ion battery in terms of transportation

Rated battery capacity: 4,5 Ah
Watt-hour rating: 162 Wh

The lithium-ion battery of the via GO V24 meet the requirements of the tests in the "UN Manual of Tests and Criteria" Part III subsection 38.3.

We hereby certify that the lithium ion battery of the via GO V24 meets the requirements of the transportation as **dangerous goods class 9**

UN 3481

under the regulations of **ADR, RID, IMDG and IATA.**

IATA UN 3481, Packing Instruction 967 PI Section I

This certificate does not apply to damaged or defect batteries according Special Provision A154, these are forbidden for transport.

Due to frequent amendments of the respective regulations this document is valid until **31st December 2023.**

Albstadt, 23rd of May 2023
Alber GmbH

A handwritten signature in blue ink, appearing to read "J. Beltman", with a long horizontal stroke extending to the right.

Johannes Beltman
Managing Director