

## Certificate in relation to air transportation by passengers with walking disabilities

We hereby certify that our product:

### Invacare Aviva FX MPS Maxx

including the non-spillable valve regulated lead acid batteries (VRLA) fulfilling special provision A67 (therefore not restricted according to UN 2800) are appropriate to be carried by passengers according

### IATA chapter 2.3.2.2

and may therefore be carried as checked baggage or in the cargo hold without restriction. Condition precedent to meet the above listed requirement is the transportation in a way that an unintentionally switching on and damage by other luggage is avoided.

Before transporting your mobility device, make sure the motors are engaged and that the remote is switched off. Invacare strongly recommends that you additionally disconnect or remove the batteries.

The pilot-in-command must be informed about the location of the battery.

We strongly recommend obtaining the approval by corresponding Airline Company before commencement date respectively departure date.

This certificate does not apply to damaged or defect batteries.

Due to the frequent amendments of the respective regulations, which partly occur on an annual basis this certificate is valid until **31<sup>st</sup> December 2026**.

Isny, 14<sup>th</sup> of January 2026

Invacare GmbH



Johannes Beltman  
Managing Director